


NEEDLE INSERTION TRAINING MODEL

Digital Twins of Medical Devices
4CBLW00-26



Group 5

1558498 Yerin Choi
1930834 Megan van Gerwen
1972529 Jelle Hiddema
2048345 Jori Schlangen
1832360 Shehab Awany
1942263 Aleksandr Vardanian

THE PROJECT

PROCEDURE

- Patient's arm is comfortably positioned, and an appropriate vein is selected.
- Needle is inserted at a shallow angle approximately 15° – 30°
- Successful insertion is indicated by a slight resistance followed by a distinct sensation "pop" or visual blood return.
- Once confirmed, blood is drawn or an IV line is established, followed by careful needle withdrawal.

POSSIBLE COMPLICATIONS

- Bruising or hematoma
- Pain or discomfort during insertion
- Difficulty in accessing the vein (multiple attempts needed)
- Patient anxiety or faintness



THE PROJECT

PROBLEM STATEMENT

- Limited hands-on practice for needle insertion
- Insufficient real-time visual feedback
- Existing simulators are expensive and not customisable

PROJECT GOAL

Develop a low-cost, customizable digital twin providing a realistic tactile experience and immediate, sensor-driven feedback.

STAKEHOLDERS

- Nurses (in training)
- Doctors (in training)

USE

Practise needle insertion technique into the lower arm (for IV and blood collection).



REQUIREMENTS

Must have

- R1.** Real-time tracking of needle pose
- R2.** Measurement of insertion force
- R3.** Physics-based tissue simulation
- R4.** User safety and reset capability


Should have

- R5.** Visual feedback interface
- R6.** Realistic physical prototype

Could have


- R7.** Exact vein location
- R8.** Precise and accurate sensors

Will not have

- R9.** Complicated needle insertion model
 - R10.** Identical syringe
- 


PHYSICAL MODEL

Final Setup




PHYSICAL MODEL

Arduino
Nano




PHYSICAL MODEL

Orientation, location and pressure




Orientation

- BNO055 motion sensor outputs Euler angles and quaternions
- Mounted on the barrel
- Visualizes needle tilt and rotation in real time
- Assists maintaining correct insertion angle (target 15° – 30°)



PHYSICAL MODEL


Orientation, location and pressure



Location

- VL6180X ToF sensor (15–200 mm range)
- Mounted under syringe barrel
- Measures distance to skin
- Needle insertion depth
- Find needle location by using orientation

TIME OF FLIGHT SENSOR



PHYSICAL MODEL

Orientation, location and pressure

Pressure


- FSR (RP-S40-ST) pressure sensor
- Detects resistance change, not exact Newton values
- Senses force changes during insertion
 - Initial contact => small force
 - Insertion => increasing force
 - Vessel entry => sudden drop
- Placed under the dummy arm



Data integration

- All sensor values are transmitted via Arduino Nano Every
 - Sent to Python-based simulation using serial communication


FINAL MODEL



SIMULATION

Deformable tissue

- Using a custom FEM (finite element method) triangular mesh
- Implemented with Young's modulus and Poisson ratio




Goal


- Adding more realism
- Improving the quality for medical training purposes

Performance issues

- Real-time deformation in the simulation required much performance
- A lot of lag when active deformation took place (dropping under 5 fps)
- Not feasible for our project's purpose and goal
- Decided to remove the component for final product



SIMULATION




SIMULATION

Color feedback



- Offers direct feedback about the needle insertion to the user of the model
- Based on the pressure sensor the tissue block turns red in the simulation
- Simulation offer real-time direct feedback

Acheived goals

- This version of the simulation runs very smooth, making it suitable for training
- By implementing feedback by color we still managed to have an effective training model with direct feedback to the user



REAL SCENARIO



LIMITATIONS

Sensors

- Fixed location insertion point
- Attachment wires on needle
- Only measure insertion depth
- Pressure measured under entire arm

Software

- Detailed OR Fast
- Still runs in SoFa
- Optimising visual feedback
- Deformable tissue had performance issues

FUTURE DEVELOPMENT

In the future, this model can be further developed to train medical students for more **complicated needle insertion procedures** on the lower arm.

The twin could also be improved by using more advanced sensors, that can track 3D location of the needle.



THE END
THANKS FOR LISTENING